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Model : NicoLet iS10 (Thermo)

I. Purpose:

1. Research on molecular structure.
2. Identification of compounds and functional groups.

II. Features of the equipment:

1. Can be measured 4000 ~ 650 cm^{-1} IR region of the spectrum. (high-sensitivity MCT Detector).
2. With ATR attachment (Attenuated Total Reflectance), the spectrum of sample surface can be measured.
3. With removable liquid sample cells equipment with temperature control device.
4. Polarizing film (manual KRS-5 Polarizer).

III. Sample preparation instructions and precautions:

1. Solid sample mixed with KBr powder in an appropriate mixing ratio (e.g.1:100) pressed into a tablet and sandwiched between two salt plates prior to the measurement by FTIR.
2. A drop of low volatility of liquid samples can be sandwiched between two pieces of salt plates (NaCl), similar to the preparation of solid samples. For samples of aqueous solution without deliquescence of non-hygroscopic window materials, such as AgCl.should be used. Do not use salt plates, such as NaCl and KBr.
3. For highly volatile liquid samples, liquid sampling cell is available to be used.
4. The ATR (Attenuated Total Reflectance) measuring attachment can be used for the analysis of the material surface.

IV. Restrictions:

1. Background spectrum of the water should be kept below 10% in order to get satisfactory FTIR spectrum.
2. Frequent entry and exit of the FTIR room should be avoided so as to avoid excessive CO_2 and H_2O contents from the atmosphere adversely influencing the FTIR spectrum of the sample.

V. Application:

1. Registration for application in advance is required. Users are required to operate the instrument themselves.
2. For users outside the department, application should be signed by the chairman of the Chemical Engineering Department and the professor-in-charge for the FTIR.

VI. Professor-In-Charge: Yan-Jyi Huang (黃延吉)